

Prism Spectrometer

Rahmanyaz Annyyev, Hikmat Gulaliyev

02 March 2024

Abstract

In this experiment, we utilize a prism spectrometer—an optical device that separates light into its constituent frequencies—to measure the refractive index of a prism. A collimator with a slit is used to produce a parallel beam of light emitted by a mercury lamp, which is then passed through the prism. The light is then refracted and dispersed into its constituent frequencies. The angle of deviation of the light is measured and used to calculate the refractive index of the prism using a special relation. A graph of the index of refraction versus wavelength is also plotted. The results are compared to the theoretical values and the limitations of the experiment are discussed.

1 Introduction

The visible spectrum is the portion of electromagnetic radiation that is visible to the human eye. It is composed of light with wavelengths between 380 and 780 nm. Each wavelength, or frequency of light is perceived as a different color by the human eye[1].

A prism spectrometer is an optical device used in spectrum analysis. It is composed of a collimator, a prism, and a telescope.

The collimator is a device that produces a parallel beam of light. It is composed of an adjustable slit on one end and a lens on the other. The slit is used to control the width of the beam of light, and the lens is used to focus the light into a parallel beam.

The prism is used to refract and disperse the light into its constituent frequencies. The angle of deviation of the light is measured and used to calculate the refractive index of the prism using a special relation.

Snell's law is given by

$$n_1 \sin(\theta_1) = n_2 \sin(\theta_2), \tag{1}$$

where n_1 and n_2 are the refractive indices of the two media, and θ_1 and θ_2 are the angles of incidence and refraction, respectively

2 Data & Results

Calculations for table 1:

Calculations for table 2:

$$\delta_{min}(\text{Violet}) = 360.00^\circ - 351.80^\circ + 29.82^\circ = 38.02^\circ \quad (2)$$

$$\delta_{min}(\text{Blue}) = 360.00^\circ - 351.80^\circ + 29.51^\circ = 38.31^\circ \quad (3)$$

$$\delta_{min}(\text{Green}) = 360.00^\circ - 351.80^\circ + 27.63^\circ = 40.19^\circ \quad (4)$$

$$\delta_{min}(\text{Orange}) = 360.00^\circ - 351.80^\circ + 28.64^\circ = 39.18^\circ \quad (5)$$

And the refractive index is given by

$$n = \frac{\sin\left(\frac{\delta_{min} + \alpha}{2}\right)}{\sin\left(\frac{\alpha}{2}\right)}, \quad (6)$$

where α is the angle of the prism. Hence, the refractive index for each wavelength is given by

$$n_{\text{Violet}} = \frac{\sin\left(\frac{38.02^\circ + 59.87^\circ}{2}\right)}{\sin\left(\frac{59.87^\circ}{2}\right)} = 1.77 \quad (7)$$

$$n_{\text{Blue}} = \frac{\sin\left(\frac{38.31^\circ + 59.87^\circ}{2}\right)}{\sin\left(\frac{59.87^\circ}{2}\right)} = 1.78 \quad (8)$$

$$n_{\text{Green}} = \frac{\sin\left(\frac{40.19^\circ + 59.87^\circ}{2}\right)}{\sin\left(\frac{59.87^\circ}{2}\right)} = 1.79 \quad (9)$$

$$n_{\text{Orange}} = \frac{\sin\left(\frac{39.18^\circ + 59.87^\circ}{2}\right)}{\sin\left(\frac{59.87^\circ}{2}\right)} = 1.78 \quad (10)$$

The angle of prism is given by

$$\alpha = \frac{\text{position 1} + \text{position 2}}{2} \quad (11)$$

Position #1	Position #2	α_i
$52^\circ 46'$	$293^\circ 01'$	59.86°
$52^\circ 52'$	$293^\circ 05'$	59.88°
$52^\circ 48'$	$293^\circ 03'$	59.86°
		$\bar{\alpha} = 59.87^\circ$

Table 1: Data for the prism angle, α .

λ (nm)	Position #1	Position #2	δ_{min}	n
404.6 Violet	351°48'	29°49'	38.02	1.77
435.8 Blue	351°48'	29°28'	38.31	1.78
546.1 Green	351°48'	27°40'	40.19	1.79
579.1 Orange	351°48'	28°41'	39.18	1.78

Table 2: Data for the refractive index versus wavelength.

3 Discussion & Conclusion

References

- [1] Robert T. Marcus. *The Measurement of Color*. Vol. 1. AZimuth. North-Holland, 1998, pp. 31–96. ISBN: 978-0-444-82825-6. DOI: [https://doi.org/10.1016/S1387-6783\(98\)80005-6](https://doi.org/10.1016/S1387-6783(98)80005-6).